



Self International System Srl
Via Vigna del Pigno, 44 - 00060 Riano - Rome - Italy

Numero Verde
800-271525

info@capribeautyline.com
www.capribeautyline.com



F
A
C
E

B
O
D
Y



ULTRA POWER TECH

SPECIFIC LINE FOR ESTHETIC DEVICE

Capri Beauty Line has created a **new high-quality line** born from the **need of our beauticians to have specific products to use with device** located in the Beauty Institutes as:

- Ultrasound Gel;
- Radio-frequency creams:
- Phialoids concentrated activating for electroporators to be added or to cream or conductive gel to enhance the desired effect.

Thanks to specially designed formulations, the products of the Ultra Power Tech Line can be transmitted through:

- ELECTROPORATORS
- ULTRASOUND
- IONTOPHORESIS
- HYPERBARIC OXYGEN

They are **dermo cosmetics** innovative and effective, to combat the major imperfections of the face and body. They contain natural active, without the use of derivatives of animals or Parabens and paraffins. All products of the line are characterized by a **high penetration index**, act in depth, and they have the right percentage of functional substances for the optimization of aesthetic results.



CONDUCTIVE CREAM FOR ESTHETIC DEVICE

COD. 57001_ML 500

Moisturizing cream particularly rich of minerals that promote its electro-conductive properties. Its suitable for energy-transfer therapies (resistive and capacitive) to be made with the appropriate devices. The product constitutes an indispensable aid for the making of any instrumental treatment and / or diathermy manual. Adding to the conductive cream Phialoids super activating anti-wrinkle, toning, draining or anti-cellulite, will be able to maximize the assets penetration with absolutely stunning results. The high conductive properties make it absolutely necessary for the use and proper functioning of the devices for beauty treatments, promoting and improving the efficiency and effectiveness of performance. Follow the manual for the application tool.



CONDUCTIVE GEL FOR ESTHETIC DEVICE

COD. 57002_ML 500

Conductive Gel Ultra Power Tech is a conductive gel carrier of active ingredients and useful for scrolling handpiece aesthetic devices. Mixed with phialoids Concentrate activators, it is suitable for various applications including slimming treatments, revitalizing, invigorating, energizing, anti-wrinkle and moisturizing. It 'a product free of Parabens and paraffins with a safe and effective formulation which, applied according to suitable treatment protocols, provides excellent results and long lasting. Follow the manual for the application tool.



DRAINING / SLIMMING CONDUCTIVE FLUID AMPOULES

COD. 57021_12 AMPOULES ML 15

Caffeine, Carnitine, Centella Asiatica, Birch, Ivy and Pineapple: a blend of natural active ingredients that, used in combination with devices and equipment, medical / aesthetic, are able to enhance law enforcement efforts and the accumulation of liquids formation of localized fat deposits, with truly amazing results. Combined with creams or conductive gels through their activities, help to reduce the swelling due to stagnation of liquids. Even with fragile capillaries allows a significant improvement of the condition and skin blemishes caused by fluid retention and cellulite.



ANTI-WRINKLE / FIRMING CONDUCTIVE FLUID AMPOULES

COD. 57025_12 AMPOULES ML 15

Ampoules concentrate formulated for use with professional devices medical / aesthetic. They contain high percentages of active principles like Hyaluronic Acid, Vitamin PP, Vitamin A / C / E, Ginseng Extract, conveyed by the device penetrate deep to the epidermis allowing to provide all nutrients and moisturizers necessary to restore the best conditions physiological.



POST TREATMENT LENITIVE SOLUTION

COD. 57011_ML 500

It has long been looking for a totally cosmetic solution to regain the balance of epidermal tissue immediately after treatment with aesthetic devices like laser, radio frequency, Electric Arc etc etc. One of the main problems of Beauty Institutes has always been to minimize the recovery period of the treated skin. The solution Post Treatment Lenitive Ultra Power Tech is able to reduce this without risk, safely and without contraindications. The active that were chosen help to calm the states of skin redness of a different nature and counteract the local inflammatory processes. The Hypericum extract, Chamomile and Aloe Vera play an 'anti-inflammatory action and contribute to the regeneration and healing of damaged tissues ensuring a gradual and long-lasting benefit. Thanks to its special formulation, it is also indicated as a nutrient for dry and stressed skins.